Indication

BeGraft

The BeGraft Stent Graft System is indicated for permanent intraluminal placement in iliac and renal arteries* for:


- restoring and improving the patency impaired by arteriosceloris or other stenoses that may be treated by stenting.
- sealing aneurysms, acute perforations, acute ruptures and fistulas.

Furthermore, the BeGraft Stent Graft System is indicated for use

- as a bridging stent graft in patients requiring endovascular abdominal aortic aneurysm repair with a fenestrated aortic endovascular graft to keep the visceral arteries perfused**.

* Sizes approved when used for permanent intraluminal placement in iliac and renal arteries: Ø 5 mm and 6 mm in lengths of 18 mm, 22 mm, 28 mm, 38 mm, and 58 mm, Ø 7 mm in lengths of 18 mm, 23 mm, 27 mm, 37 mm, and 57 mm, Ø 8 mm - 10 mm in lengths of 27 mm, 37 mm, and 57 mm

**Sizes approved when used as a bridging stent graft for FEVAR: Ø 5 mm and 6 mm in lengths of 22 mm, 28 mm, and 38 mm, Ø 7 mm in lengths of 23 mm, 27 mm, and 37 mm, Ø 8 mm - 10 mm in lengths of 27 mm and 37 mm



Bentley

Product Specification

BeGraft



Expanded Stent Graft Diameter	Nominal Stent Graft Length (mm)	Introducer Sheath Compatibility	Catheter Lengths
5mm	18, 22*, 28*, 38*, 58	6F	75cm/ 120cm
6mm	18, 22*, 28*, 38*, 58	6F	75cm/ 120cm
7mm	18, 23*, 27*, 37*, 57	6F	75cm/ 120cm
8mm	27*, 37*, 57	6F	75cm/ 120cm
9mm	27*, 37*, 57	7F	75cm/ 120cm
10mm	27*, 37*, 57	7F	75cm/ 120cm

* Sizes indicated for FEVAR indication

Bentley

Product Specification – Iliac and Renal Indication

BeGraft



Expanded Stent Graft Diameter	Nominal Stent Graft Length (mm)	Introducer Sheath Compatibility	Catheter Lengths
5mm	18, 22, 28, 38, 58	6F	75cm/ 120cm
6mm	18, 22, 28, 38, 58	6F	75cm/ 120cm
7mm	18, 23, 27, 37, 57	6F	75cm/ 120cm
8mm	27, 37, 57	6F	75cm/ 120cm
9mm	27, 37, 57	7F	75cm/ 120cm
10mm	27, 37, 57	7F	75cm/ 120cm

Bentley

Product Specification – FEVAR Indication

BeGraft




Expanded Stent Graft Diameter	Nominal Stent Graft Length (mm)	Introducer Sheath Compatibility	Catheter Lengths
5mm	22, 28, 38	6F	75cm/ 120cm
6mm	22, 28, 38	6F	75cm/ 120cm
7mm	23, 27, 37	6F	75cm/ 120cm
8mm	27, 37	6F	75cm/ 120cm
9mm	27, 37	7F	75cm/ 120cm
10mm	27, 37	7F	75cm/ 120cm

Bentley

Product Specification

BeGraft

Stent Graft




Bentley

Product Specification

BeGraft

Balloon Catheter



Bentley